

Filter P1D09




- Conservative therapy
- Prosthodontics
- Endodontics
- Periodontics
- Implantology
- Jaw surgery
- Aesthetic surgery

Laser and application

Typical laser	Wavelength	OD	CE protection rating
Er:YAG	2940nm	4+	DI LB4+R LB3Y+M LB3
CO ₂ laser	10,600nm	6+	DI LB4+R LB3Y+M LB4Y

Filter properties

Filter function: Full protection
Filter colour: Grey-green
Visual brightness: Very good
VLT (approx.): 50%
Colour view: Very good



Available laser protection eyewear

F18
Page 14  Art. number: **F18.P1D09**

F22
Page 15  Art. number: **F22.P1D09**

F42
Page 15  Art. number: **F42.P1D09**

F27
Page 16  Art. number: **F27.P1D09**

F46
Page 18  Art. number: **F46.P1D09**



Filter P1P10




- Conservative therapy
- Prosthodontics
- Endodontics
- Periodontics
- Implantology
- Jaw surgery
- Aesthetic surgery

Laser and application

Typical laser	Wavelength	OD	CE protection rating
Nd:YAG	1064nm	8+	DLB6+IRM LB8
CO ₂ laser	10,600nm	4+	DI LB4+RLB3Y+M LB5

Filter properties

Filter function: Full/alignment protection
Filter colour: Blue
Visual brightness: Sufficient
VLT (approx.): 16%
Colour view: Good



Available laser protection eyewear

F18
Page 14  Art. number: **F18.P1P10**

F22
Page 15  Art. number: **F22.P1P10**

F42
Page 15  Art. number: **F42.P1P10**

F27
Page 16  Art. number: **F27.P1P10**

F29
Page 21  Art. number: **F29.P1P10**

F46
Page 18  Art. number: **F46.P1P10**

F47
Page 22  Art. number: **F47.P1P10**

