



Filter P1P18




- Endodontics
- Periodontics
- Soft tissue surgery
- Therapies
- Biostimulation and whitening

Laser and application

Typical laser	Wavelength	OD	CE protection rating
Diode	405-460nm	8+	D LB6 + IR LB8 + M LB8Y
Diode	810nm	6+	DIRM LB6
Diode	940+980nm	10+	D LB6 + IR LB8 + M LB8Y
Nd:YAG	1064nm	10+	D LB6 + IR LB8 + M LB10

Filter properties

Filter function: Full protection
Filter colour: Gold
Visual brightness: Very good
VLT (approx.): 47%
Colour view: Good



Available laser protection eyewear

F18 Page 14		Art. number: F18.P1P18
F22 Page 15		Art. number: F22.P1P18
F42 Page 15		Art. number: F42.P1aP18
F27 Page 16		Art. number: F27.P1P18
F29 Page 21		Art. number: F29.P1P18
F46 Page 18		Art. number: F46.P1P18



Filter P1M03




- Conservative therapy
- Prosthodontics
- Endodontics
- Periodontics
- Implantology
- Jaw surgery
- Aesthetic surgery

Laser and application

Typical laser	Wavelength	OD	CE protection rating
Nd:YAG	1064nm	8+	DLB6+IRLB8+MLB8Y
Er:YAG	2940nm	4+	DILB4+RLB3Y+MLB4

Filter properties

Filter function: Full protection
Filter colour: Grey
Visual brightness: Very good
VLT (approx.): 25%
Colour view: Very good



Available laser protection eyewear

F18 Page 14		Art. number: F18.P1M03
F22 Page 15		Art. number: F22.P1M03
F42 Page 15		Art. number: F42.P1M03
F27 Page 16		Art. number: F27.P1M03
F29 Page 21		Art. number: F29.P1M03
F46 Page 18		Art. number: F46.P1M03
F47 Page 22		Art. number: F47.P1M03

