



Filter P1M03




- Removal/treatment of skin lesions
- Removal/treatment of epidermal lesions

Laser and application

Typical laser	Wavelength	OD	CE protection rating
Nd:YAG	1064nm	8+	DLB6+IRLB8+MLB8Y
Er:YAG	2940nm	4+	DI LB4 + R LB3Y + M LB4

Filter properties

Filter function: Full protection
Filter colour: Grey
Visual brightness: Very good
VLT (approx.): 25%
Colour view: Very good



Available laser protection eyewear

F18
Page 14
 Art. number: **F18.P1M03**

F22
Page 15
 Art. number: **F22.P1M03**

F42
Page 15
 Art. number: **F42.P1M03**

F27
Page 16
 Art. number: **F27.P1M03**

F29
Page 21
 Art. number: **F29.P1M03**

F46
Page 18
 Art. number: **F46.P1M03**

F47
Page 22
 Art. number: **F47.P1M03**

For further information and to view the technical specifications of these filters, please scan the QR codes provided.



Filter P1D09




- Removal/treatment of skin lesions
- Removal/treatment of epidermal lesions

Laser and application

Typical laser	Wavelength	OD	CE protection rating
Er:YAG	2940nm	4+	DI LB4 + R LB3Y + M LB3

Filter properties

Filter function: Full protection
Filter colour: Grey-green
Visual brightness: Very good
VLT (approx.): 50%
Colour view: Very good



Available laser protection eyewear

F18
Page 14
 Art. number: **F18.P1D09**

F22
Page 15
 Art. number: **F22.P1D09**

F42
Page 15
 Art. number: **F42.P1D09**

F27
Page 16
 Art. number: **F27.P1D09**

F46
Page 18
 Art. number: **F46.P1D09**

